


# Best Seller

## 2~8°C Pharmacy Refrigerator

Intelligent Control & Constant Temperature  
Multiple Alarms Bring More Security



**YC-130L**



Temp. Uniformity



Pulling Down to 5°C



Temp. Recovery After Door Opened for 1 Min

*\*Typical data based on internal prototype testing with set point at 5°C and ambient temperature at 25°C, may vary between units.*

## Features

**Temperature Control System**  
Display temp. in 0.1°C increments

**Compact and Fit**  
It fits under standard 34" high ADA compliant countertop


**HC Nature Hydrocarbon**  
Environment-friendly Natural Gas



**Backup Battery**  
Can maintain real-time temperature display and alarms in case of power failure

**Double-pane Glass Door**  
Automatic electric heating function

**Data Logger**  
It can log 100000pcs data, 12 months of data can be exported via USB port

## Structure





| 2~8°C Pharmacy Refrigerator   |   |               |                  |
|-------------------------------|---|---------------|------------------|
| Model                         | YC-130L   |               |                  |
| Capacity                      | 130L/4.59Cu.Ft  |               |                  |
| Internal Size (W*D*H) mm/Inch | 554*510*588/21.81*20.08*23.15                           |               |                  |
| External Size (W*D*H) mm/Inch | 650*625*810/25.59*24.60*23.15                           |               |                  |
| Package Size (W*D*H) mm/Inch  | 720*703*880/28.35*27.68*34.65                           |               |                  |
| NW(Kgs/Lbs)                   | 51/112.44   |               |                  |
| GW(Kgs/Lbs)                   | 61/141.1  |               |                  |
| Performance                   |   |               |                  |
| Temperature Range             | 2~8°C   |               |                  |
| Ambient Temperature           | 16~32°C   |               |                  |
| Cooling Performance           | 5°C   |               |                  |
| Climate Class                 | N   |               |                  |
| Controller                    | Microprocessor  |               |                  |
| Display                       | Digital Display   |               |                  |
| Refrigeration                 |   |               |                  |
| Compressor                    | 1pc   |               |                  |
| Cooling Method                | Air Cooling   |               |                  |
| Defrost Mode                  | Automatical   |               |                  |
| Refrigerant                   | R290  | R600a         | R134a            |
| Insulation Thickness(mm)      | D:55, L/R:48, B:50                                      |               |                  |
| Construction                  |   |               |                  |
| External Material             | PCM   |               |                  |
| Inner Material                | Sprayed Aluminum Plate                                  |               |                  |
| Shelves                       | 3   |               |                  |
| Access Test Port              | 1pc. Ø23mm  |               |                  |
| Casters                       | 2+(2 Leveling Feet)                                     |               |                  |
| Data Logging/Export           | Record every 10 mins/Max.12 months via USB port         |               |                  |
| Interface                     | Remote Alarm Contact, RS485                             |               |                  |
| Backup Battery                | Yes   |               |                  |
| Alarms                        |   |               |                  |
| Temperature                   | High/Low Temperature, High Ambient Temperature          |               |                  |
| Electrical                    | Power Failure, Low Battery, Sensor Error                |               |                  |
| System                        | Communication error,Door ajar,Data logger not connected |               |                  |
| Electrical                    |   |               |                  |
| Power Supply(V/Hz)            | 110V/60Hz   | 220~240V;50Hz | 220~240V;50/60Hz |
| Rated Current(A)              | 2.28A   | 0.98A         | 1.1A             |

\* Zhongke Meiling reserves the right to change products and specifications without prior notice.